

PACO33952

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

Participates in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyzes dithiol-disulfide exchange reactions.

**Gene ID:**

trxA

**Uniprot**

P0AA25

**Synonyms:**

Thioredoxin 1 (Trx-1), trxA, fipA tsnC

**Immunogen:**


Recombinant Escherichia coli Thioredoxin-1 protein (2-109AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: trxA antibody at 2 $\mu$ g/ml + DH5alpha whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 12 kDa

Observed band size: 12 kDa

.