

MT2731 Antibody



PACO33976

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Antitoxin component of a toxin-antitoxin (TA) module. Neutralizes the effect of cognate toxin MT2730.

Gene ID:

MT2731

Uniprot

P9WJ10

Synonyms:

Antitoxin MT2731

Immunogen:

Recombinant Mycobacterium tuberculosis Antitoxin MT2731 protein (1-81AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A