

# ydhF Antibody



PACO33980

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May function as oxidoreductase.

**Gene ID:**

ydhF

**Uniprot**

P76187

**Synonyms:**

Oxidoreductase YdhF (EC 1. -. -. -), ydhF

**Immunogen:**

Recombinant Escherichia coli Oxidoreductase YdhF protein (1-298AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A