

PACO33996

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the biosynthesis of alanine.

Gene ID:

alaC

Uniprot

P77434

Synonyms:

Glutamate-pyruvate aminotransferase AlaC (EC 2.6.1.2), alaC, yfdZ

Immunogen:

Recombinant Escherichia coli Glutamate-pyruvate aminotransferase AlaC protein (1-412AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A