

ALTA1 Antibody



PACO34020

Product Information

Size:

50ug

Reactivity:

Alternaria alternata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ALTA1

Uniprot

P79085

Synonyms:

Major allergen Alt a 1 (allergen Alt a 1), ALTA1

Immunogen:

Recombinant Alternaria alternata Major allergen Alt a 1 protein (19-157AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A