

Defensin D2 Antibody



PACO34028

Product Information

Size:

50ug

Reactivity:

Spinacia oleracea

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Antimicrobial peptide. Active against Fusarium spp. , Gram-positive and Gram-negative bacterial pathogens.

Gene ID:**Uniprot**

P81571

Synonyms:

Defensin D2 (Antimicrobial peptide D2) (So-D2) (Fragment)

Immunogen:

Recombinant Spinacia oleracea Defensin D2 protein (1-52AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A