

PACO34074

Product Information

Size:

50ug

Reactivity:

Pseudomonas aeruginosa

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

D-galactose specific lectin. Binds in decreasing order of affinity: melibiose, methyl-alpha-D-galactoside, D-galactose, methyl-beta-D-galactoside, N-acetyl-D-galactosamine. Similar to plant lectins in its selective (carbohydrate-specific) hemagglutinating activity.

Gene ID:

lecA

Uniprot

Q05097

Synonyms:

PA-I galactophilic lectin (PA-IL) (Galactose-binding lectin), lecA, pa1L

Immunogen:

Recombinant Pseudomonas aeruginosa PA-I galactophilic lectin protein (2-122AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A