

PACO34090

Product Information

Size:

50ug

Reactivity:

Zaire ebolavirus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GP

Uniprot

P0C773

Synonyms:

Super small secreted glycoprotein (SsGP), GP

Immunogen:

Recombinant Zaire ebolavirus Super small secreted glycoprotein (33-298AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A