

# p46 Antibody



PACO34094

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mycoplasma hyopneumoniae

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

p46

**Uniprot**

P0C0J8

**Synonyms:**

46 kDa surface antigen (p46), p46

**Immunogen:**

Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A