

p46 Antibody



PACO34098

Product Information

Size:

50ug

Reactivity:

Mycoplasma hyopneumoniae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

p46

Uniprot

P0C0J8

Synonyms:

46 kDa surface antigen (p46), p46

Immunogen:

Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A