

# IktA Antibody



PACO34106

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mannheimia haemolytica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IktA

**Uniprot**

P0C085

**Synonyms:**

Leukotoxin (Lkt), IktA

**Immunogen:**

Recombinant Mannheimia haemolytica Leukotoxin protein (715-953AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A