

PACO34110

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Efficient electron donor for the essential enzyme ribonucleotide reductase. Is also able to reduce the interchain disulfide bridges of insulin.

**Gene ID:**

trxC

**Uniprot**

P0AGG4

**Synonyms:**

Thioredoxin 2 (Trx-2) (EC 1.8.1.8) (Protein-disulfide reductase), trxC, yfiG

**Immunogen:**

Recombinant Escherichia coli Thioredoxin 2 protein (1-139AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A