

# PMP2 Antibody



PACO34114

---

## Product Information

**Size:**

50ug

**Reactivity:**

Horse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in lipid transport protein in Schwann cells. May bind cholesterol.

**Gene ID:**

PMP2

**Uniprot**

P0C6G6

**Synonyms:**

Myelin P2 protein, PMP2

**Immunogen:**

Recombinant Horse Myelin P2 protein (2-132AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A