

PMP2 Antibody, FITC conjugated



PACO34116

Product Information

Size:

50ug

Reactivity:

Horse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in lipid transport protein in Schwann cells. May bind cholesterol.

Gene ID:

PMP2

Uniprot

P0C6G6

Synonyms:

Myelin P2 protein, PMP2

Immunogen:

Recombinant Horse Myelin P2 protein (2-132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A