

PACO34118

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

The biological activity of the toxin is produced by the A chain, which activates intracellular adenyl cyclase.

Gene ID:

eltB

Uniprot

P0CK94

Synonyms:

Heat-labile enterotoxin B chain (LT-B, human) (LTH-B), eltB, ltpB

Immunogen:

Recombinant Escherichia coli Heat-labile enterotoxin B chain protein (22-124AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A