

PACO34122

Product Information

Size:

50ug

Reactivity:

Gallus gallus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Seems to have at least five active inhibitory domains; two for trypsin, two for chymotrypsin and one for elastase.

Gene ID:

OIH

Uniprot

P10184

Synonyms:

Ovoinhibitor, OIH

Immunogen:

Recombinant Gallus gallus Ovoinhibitor protein (26-472AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A