

# Mup2 Antibody



PACO34158

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds pheromones that are released from drying urine of males. These pheromones affect the sexual behavior of females.

**Gene ID:**

Mup2

**Uniprot**

P11589

**Synonyms:**

Major urinary protein 2 (MUP 2), Mup2

**Immunogen:**

Recombinant Mouse Major urinary protein 2 protein (19-180AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A