

GRA1 Antibody



PACO34166

Product Information

Size:

50ug

Reactivity:

Toxoplasma gondii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GRA1

Uniprot

P13403

Synonyms:

Dense granule protein 1 (Protein GRA 1) (Major antigen p24), GRA1

Immunogen:

Recombinant Toxoplasma gondii Dense granule protein 1 protein (26-187AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A