

PACO34182

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Conversion of glycerol 3-phosphate to dihydroxyacetone. Uses molecular oxygen or nitrate as electron acceptor.

**Gene ID:**

glpD

**Uniprot**

P13035

**Synonyms:**

Aerobic glycerol-3-phosphate dehydrogenase (EC 1.1.5.3), glpD, glyD

**Immunogen:**

Recombinant Escherichia coli Aerobic glycerol-3-phosphate dehydrogenase protein (1-501AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A