

AP24 Antibody



PACO34190

Product Information

Size:

50ug

Reactivity:

Nicotiana tabacum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

AP24

Uniprot

P14170

Synonyms:

Osmotin, AP24


Immunogen:

Recombinant Nicotiana tabacum Osmotin protein (22-246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: *Nicotiana tabacum* (40 µg; g, 20 µg; g, 10 µg; g). All lanes: AP24 antibody at 3 µg/ml. Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 27 kDa. Observed band size: 27 kDa.