

PACO34214

Product Information

Size:

50ug

Reactivity:

Bordetella pertussis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Agglutinin that binds to eukaryotic cells; a process mediated by the R-G-D sequence. Pertactin may have a role in bacterial adhesion, and thus play a role in virulence. May contribute to the disease state of whooping cough.

Gene ID:

prn

Uniprot

P14283

Synonyms:

Pertactin autotransporter (P.93) [Cleaved into: Outer membrane protein P.69; Pertactin translocator], prn, omp69A

Immunogen:

Recombinant Bordetella pertussis Pertactin autotransporter protein (632-910AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A