

PACO34218

Product Information

Size:

50ug

Reactivity:

Epstein-Barr virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SCP

Uniprot

P14348

Synonyms:

Small capsomere-interacting protein, SCP

Immunogen:

Recombinant Epstein-Barr virus Small capsomere-interacting protein (31-176AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A