

PACO34234

Product Information

Size:

50ug

Reactivity:

Streptomyces hygroscopicus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This enzyme is an effector of phosphinothricin tripeptide (PTT or bialaphos) resistance. Inactivates PTT by transfer of an acetyl group.

Gene ID:

bar

Uniprot

P16426

Synonyms:

Phosphinothricin N-acetyltransferase (PPT N-acetyltransferase) (EC 2.3.1.183)
(Phosphinothricin-resistance protein), bar

Immunogen:

Recombinant Streptomyces hygroscopicus Phosphinothricin N-acetyltransferase protein (1-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A