

toxA Antibody



PACO34290

Product Information

Size:

50ug

Reactivity:

Peptoclostridium difficile

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Only after the enteral delivery of the enterotoxin A may the characteristic disease called pseudomembranous colitis be induced.

Gene ID:

toxA

Uniprot

P16154

Synonyms:

Toxin A (EC 3.4.22. -), toxA, tcdA

Immunogen:

Recombinant Peptoclostridium difficile Toxin A protein (2387-2710AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A