

Prothoracicotropic hormone Antibody



PACO34298

Product Information

Size:

50ug

Reactivity:

Bombyx mori

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

PTTH is a brain secretory polypeptide of insects which stimulates the prothoracic glands to produce and release ecdysone, the steroid essential to insect development. Peptides P2K and P6K are presumed to be cleaved post-translationally and may play some unknown physiologically or developmentally important functions.

Gene ID:**Uniprot**

P17219

Synonyms:

Prothoracicotropic hormone (PTTH) [Cleaved into: P2K; P6K; Prothoracicotropic hormone]

Immunogen:

Recombinant Bombyx mori Prothoracicotropic hormone protein (116-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A