

Nuclease P1 Antibody



PACO34370

Product Information

Size:

50ug

Reactivity:

Penicillium citrinum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Hydrolyzes only single-stranded DNA and RNA without apparent specificity for bases.

Gene ID:

Uniprot

P24289

Synonyms:

Nuclease P1 (EC 3.1.30.1) (Deoxyribonuclease P1) (Endonuclease P1)

Immunogen:

Recombinant Penicillium citrinum Nuclease P1 protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A