

# cfaE Antibody



PACO34374

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

cfaE

**Uniprot**

P25734

**Synonyms:**

CFA/I fimbrial subunit E (Colonization factor antigen I subunit E), cfaE

**Immunogen:**

Recombinant Escherichia coli CFA/I fimbrial subunit E protein (1-360AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A