

Pectate lyase 5 Antibody, HRP conjugated



PACO34387

Product Information

Size:

50ug

Reactivity:

Ambrosia artemisiifolia

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has pectate lyase activity.

Gene ID:**Uniprot**

P27759

Synonyms:

Pectate lyase 5 (EC 4.2.2.2) (Antigen Amb a I) (Antigen E) (AgE) (Pollen allergen Amb a 1.1) (allergen Amb a 1.1)

Immunogen:

Recombinant Ambrosia artemisiifolia Pectate lyase 5 protein (26-396AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A