

PACO34434

---

## Product Information

**Size:**

50ug

**Reactivity:**

Aspergillus niger

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the hydrolysis of inorganic orthophosphate from phytate.

**Gene ID:**

phyA

**Uniprot**

P34752

**Synonyms:**

3-phytase A (EC 3.1.3.8) (3 phytase A) (Myo-inositol hexakisphosphate phosphohydrolase A) (Myo-inositol-hexaphosphate 3-phosphohydrolase A), phyA

**Immunogen:**

Recombinant Aspergillus niger 3-phytase A protein (24-467AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A