

phyA Antibody, Biotin conjugated



PACO34437

Product Information

Size:

50ug

Reactivity:

Aspergillus niger

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the hydrolysis of inorganic orthophosphate from phytate.

Gene ID:

phyA

Uniprot

P34752

Synonyms:

3-phytase A (EC 3.1.3.8) (3 phytase A) (Myo-inositol hexakisphosphate phosphohydrolase A) (Myo-inositol-hexaphosphate 3-phosphohydrolase A), phyA

Immunogen:

Recombinant Aspergillus niger 3-phytase A protein (24-467AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A