

gluD Antibody, Biotin conjugated



PACO34445

Product Information

Size:

50ug

Reactivity:

Clostridium difficile

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

gluD

Uniprot

P27346

Synonyms:

NAD-specific glutamate dehydrogenase (NAD-GDH) (EC 1.4.1.2), gluD

Immunogen:

Recombinant Clostridium difficile NAD-specific glutamate dehydrogenase protein (1-421AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A