

SUS1 Antibody, Biotin conjugated



PACO34449

Product Information

Size:

50ug

Reactivity:

Saccharum officinarum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Sucrose-cleaving enzyme that provides UDP-glucose and fructose for various metabolic pathways.

Gene ID:

SUS1

Uniprot

P31925

Synonyms:

Sucrose synthase (EC 2.4.1.13) (Sucrose-UDP glucosyltransferase) (Fragment), SUS1

Immunogen:

Recombinant Saccharum officinarum Sucrose synthase protein (1-218AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A