

PACO34450

Product Information

Size:

50ug

Reactivity:

Saccharomyces cerevisiae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor binding subunit of the a-agglutinin heterodimer. S. cerevisiae a and alpha cells express the complementary cell surface glycoproteins a-agglutinin and alpha-agglutinin, respectively, which interact with one another to promote cellular aggregation during mating.

Gene ID:

AGA2

Uniprot

P32781

Synonyms:

A-agglutinin-binding subunit, AGA2

Immunogen:

Recombinant Saccharomyces cerevisiae A-agglutinin-binding subunit protein (19-87AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A