

PACO34458

Product Information

Size:

50ug

Reactivity:

Entamoeba histolytica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the reversible phosphorylation of pyruvate and phosphate. In *E. histolytica* and *C. symbiosus*, PPK functions in the direction of ATP synthesis.

Gene ID:

PPDK

Uniprot

P37213

Synonyms:

Pyruvate, phosphate dikinase (EC 2.7.9.1) (Pyruvate, orthophosphate dikinase), PPK

Immunogen:

Recombinant *Entamoeba histolytica* Pyruvate, phosphate dikinase protein (1-342AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A