

Polygalacturonase Antibody



PACO34462

Product Information

Size:

50ug

Reactivity:

Cryptomeria japonica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P43212

Synonyms:

Polygalacturonase (PG) (EC 3.2.1.15) (Allergen Cry j II) (Major pollen allergen Cry j 2) (Pectinase) (allergen Cry j 2)

Immunogen:

Recombinant Cryptomeria japonica Polygalacturonase protein (46-433AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A