

AMA-1 Antibody



PACO34466

Product Information

Size:

50ug

Reactivity:

Plasmodium falciparum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in parasite invasion of erythrocytes.

Gene ID:

AMA-1

Uniprot

P22621

Synonyms:

Apical membrane antigen 1 (Merozoite surface antigen), AMA-1, PF83

Immunogen:

Recombinant Plasmodium falciparum Apical membrane antigen 1 protein (25-546AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A