

PACO34478

Product Information

Size:

50ug

Reactivity:

Porphyromonas gingivalis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thiol protease which is believed to participate in intracellular degradation and turnover of proteins. Its proteolytic activity is a major factor in both periodontal tissue destruction and in bacterial host defense mechanisms. Activates complement C3 and C5.

Gene ID:

rgpA

Uniprot

P28784

Synonyms:

Gingipain R1 (EC 3.4.22.37) (Arg-gingipain) (Gingipain 1) (RGP-1), rgpA, rgp1

Immunogen:

Recombinant Porphyromonas gingivalis Gingipain R1 protein (228-720AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A