

uspF Antibody



PACO34502

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

ATP binding, cell adhesion, response to oxidative stress.

Gene ID:

uspF

Uniprot

P37903

Synonyms:

Universal stress protein F, uspF, ynaF yzzL

Immunogen:

Recombinant Escherichia coli Universal stress protein F protein (1-144AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A