

rodA Antibody



PACO34506

Product Information

Size:

50ug

Reactivity:

Neosartorya fumigata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cell wall protein regularly arranged in interwoven fascicules of clustered proteinaceous microfibrils, or rodlets, to form the outer spore coat protein. It is involved in resistance to environmental stress and may well be associated with conidial hydrophobicity. It is important in the morphogenesis of the dispersible conidia.

Gene ID:

rodA

Uniprot

P41746

Synonyms:

Hydrophobin (Rodlet protein), rodA, hyp1

Immunogen:

Recombinant Neosartorya fumigata Hydrophobin protein (19-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A