

pal Antibody



PACO34510

Product Information

Size:

50ug

Reactivity:

Legionella pneumophila

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Very strongly associated with the peptidoglycan.

Gene ID:

pal

Uniprot

P26493

Synonyms:

Peptidoglycan-associated lipoprotein (19 kDa surface antigen) (PPL), pal, ppIA

Immunogen:

Recombinant Legionella pneumophila Peptidoglycan-associated lipoprotein protein (22-176AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A