

# BBE1 Antibody



PACO34518

---

## Product Information

**Size:**

50ug

**Reactivity:**

Eschscholzia californica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Essential to the formation of benzophenanthridine alkaloids in the response of plants to pathogenic attack. Catalyzes the stereospecific conversion of the N-methyl moiety of (S)-reticuline into the berberine bridge carbon of (S)-scoulerine.

**Gene ID:**

BBE1

**Uniprot**

P30986

**Synonyms:**

Reticuline oxidase (EC 1.21.3.3) (Berberine bridge-forming enzyme) (BBE) (Tetrahydroprotoberberine synthase), BBE1

**Immunogen:**

Recombinant Eschscholzia californica Reticuline oxidase protein (24-538AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A