

OV16 Antibody



PACO34522

Product Information

Size:

50ug

Reactivity:

Onchocerca volvulus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

OV16

Uniprot

P31729

Synonyms:

OV-16 antigen, OV16

Immunogen:

Recombinant Onchocerca volvulus OV-16 antigen protein (17-197AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A