

# MSA2 Antibody



PACO34534

---

## Product Information

**Size:**

50ug

**Reactivity:**

Plasmodium falciparum

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May play a role in the merozoite attachment to the erythrocyte.

**Gene ID:**

MSA2

**Uniprot**

P50498

**Synonyms:**

Merozoite surface antigen 2 (MSA-2) (45 kDa merozoite surface antigen), MSA2

**Immunogen:**

Recombinant Plasmodium falciparum Merozoite surface antigen 2 protein (109-246AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A