

# TES-26 Antibody



PACO34570

---

## Product Information

**Size:**

50ug

**Reactivity:**

Toxocara canis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TES-26

**Uniprot**

P54190

**Synonyms:**

26 kDa secreted antigen (Toxocara excretory-secretory antigen 26) (TES-26), TES-26

**Immunogen:**

Recombinant Toxocara canis 26 kDa secreted antigen protein (22-262AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A