

TES-26 Antibody, Biotin conjugated



PACO34573

Product Information

Size:

50ug

Reactivity:

Toxocara canis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TES-26

Uniprot

P54190

Synonyms:

26 kDa secreted antigen (Toxocara excretory-secretory antigen 26) (TES-26), TES-26

Immunogen:

Recombinant Toxocara canis 26 kDa secreted antigen protein (22-262AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A