

PACO34578

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

CMC4

**Uniprot**

P56277

**Synonyms:**

Cx9C motif-containing protein 4 (Mature T-cell proliferation 1 neighbor protein) (Mature T-cell proliferation-1 type A) (MTCP-1 type A) (Protein p8 MTCP-1) (p8MTCP1), CMC4, C6.1B MTCP1 MTCP1NB

**Immunogen:**


Recombinant Human Cx9C motif-containing protein 4 protein (1-68AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO34578 at dilution of 1:100.