

gag Antibody, HRP conjugated



PACO34643

Product Information

Size:

50ug

Reactivity:

Moloney murine leukemia virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

gag

Uniprot

P03332

Synonyms:

Gag polyprotein (Pr65gag) (Core polyprotein) [Cleaved into: Matrix protein p15 (MA); RNA-binding phosphoprotein p12 (pp12); Capsid protein p30 (CA); Nucleocapsid protein p10-Gag (NC-gag)], gag

Immunogen:

Recombinant Moloney murine leukemia virus Gag polyprotein protein (216-478AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A