

# NP Antibody



PACO34666

---

## Product Information

**Size:**

50ug

**Reactivity:**

Influenza B virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NP

**Uniprot**

P04666

**Synonyms:**

Nucleoprotein (Nucleocapsid protein) (Protein N), NP

**Immunogen:**

Recombinant Influenza B virus Nucleoprotein protein (1-168AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A