

PACO34682

Product Information

Size:

50ug

Reactivity:

Bacillus thuringiensis subsp

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes colloidosmotic lysis by binding to the midgut epithelial cells of many lepidopteran larvae.

Gene ID:

cry1Ac

Uniprot

P05068

Synonyms:

Pesticidal crystal protein Cry1Ac (133 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIA(c)), cry1Ac, cry1A(c) cry218 cryIA(c)

Immunogen:

Recombinant Bacillus thuringiensis subsp Pesticidal crystal protein cry1Ac protein (972-1178AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A