

# ompC Antibody



PACO34698

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Forms pores that allow passive diffusion of small molecules across the outer membrane.

**Gene ID:**

ompC

**Uniprot**

P06996

**Synonyms:**

Outer membrane protein C (Outer membrane protein 1B) (Porin OmpC), ompC, meoA par

**Immunogen:**

Recombinant Escherichia coli Outer membrane protein C protein (22-367AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A